



Features:-

- > High Stability and Accuracy
- > 3.5-digit LED display
- > Slope Control
- > Proper Tested



Digital pH Meter S-901/S-902

These pH meter are on the top of people's priority due to their fantastic quality and high durability. They are fabricated with the help of perfect raw materials and international standard technology. Our engineers R & D Team and experts' manufacturer digital pH Meter of complete international standards, customized solutions and customer centric approach is our top priority and business policy.

Technical Specifications : -

Model No.	S-901	S-902
Display	3.5-digit LED	4-digit LED
Resolution	0.01 pH, 1 mV	0.01 pH, 1 mV
Range	0 to 14.00 pH	0 to 14.00 pH, 0 to ± 1999 mV
Calibration	2 Point (4pH, 7pH)	3 Point (4pH, 7pH and 9.20pH)
Accuracy	± 0.0 pH, ± 1 mV	± 0.0 pH, ± 1 mV
Temperature	Manual 0 to 100 Degree Celsius	Auto 0 to 100 Degree Celsius
Slope Correction	80 to 120%	80 to 120% Auto Ranging
Power Supply	230V, $\pm 10\%$ AC, 50Hz	230V, $\pm 10\%$ AC, 50Hz
Dimension	280X180X88 mm (LXBXH)	195X225X65 mm (LXBXH)
Weight	3 KG	2 KG
Accessories	Combination pH Electrode, Buffer Tablets (4pH, 7pH), pH Electrode Stand, Instruction Manual, Dust Cover	Combination pH Electrode, Buffer Tablets (4pH, 7pH, 9.20 pH), pH Electrode Stand, Instruction Manual, Dust Cover, Temp. probe.