

This Microscope comes with Achromatic optical system, which provides a very good optical quality image. Also, with 5 Binocular head system, provided with a wooden table, makes it easier to simultaneously see the sample and study them together for 5 people. High intensity LED lighting makes it easier to view the samples without any stress. This makes the microscope optimum for Medical & Dental College Students.

- > 5 Binocular head for multiple view
- > Wide Field Eyepiece WF 10x/18mm
- > Quadruple Nosepiece
- > Achromatic Objective 4x, 10x, 40x, 100x
- > Swing Condenser N.A.0.9/0.25
- > Coaxial Coarse and Fine Adjustment,
- > Fine Division 0.002mm
- > Plain Layer Mechanical Stage 140mm x 140mm
- > Halogen lamp 6V/30W, Brightness Adjustable
- > LED Pointer, Brightness Adjustable

Packaging: Inner Thermocol/Styrofoam Packing and Outer Cardboard Box.

PENTAHEAD MICROSCOPE

X-555

