



# *ultra* Series

- > Monocular Head having Hard coated prisms
- > Quadruple revolving Ball Bearing Nose piece
- > Co-axial Focussing system
- > Double layered mechanical stage with low positioned coaxial knobs
- > Bright LED Light with Solid state power supply circuitry.
- > Plan Achromatic objectives, Extra Wide Field Eye Piece {Hard Coated} along with high quality condenser lens for better results for all magnifications.

**Objectives:** 4x, 10x, 40x, 100x (oil)

**Eye pieces :** EWF 10x (18mm Paired)

**Packaging:** Inner Thermocol/Styrofoam Packing and Outer Cardboard Box.

**Optional Packaging:**

1) Wooden Box outer cabinet with lock and key (Chargeable).



## MONOCULAR MICROSCOPE BM-8mo

*ultra*

