

Projection Microscope

(KW-500)

Illumination

It comes equipped with 12V 100W Halogen Lamp

Stage

110 x 110mm square stage is provided with stage clips to hold the slide (optional Mechanical Stage)

Focussing

Focussing control knob on both sides of the instrument is conveniently located for easy operation of the rack and pinion with 15mm range.

Nose Piece

Triple revolving nosepiece with positive click stops.

Screen & Wall Projection

160mm dia. optically true ground glass screen is provided to have high resolving power. The screen housing is hinged to the body frame and can be removed to one side to have projection on the wall or screen in a dark room

Objectives

This Microscope is equipped with three Achromatic objectives 5x, 10x & 20x

Packaging: Inner Thermocol/Styrofoam Packing and Outer Cardboard Box.

