

## DS-50

The excellent appearance  
in accord with ergonomic design

Integrated body design  
makes products maintains  
long-term stability.

High quality 3W LED lamp  
and long usage life.

- > Binocular Head Inclined at 45°, Interpupillary 55-75mm
- > Wide Field Eyepiece EW 10x & 18x
- > Objective - 2x/4x or 1x/3x
- > Working Distance 100mm
- > Stage 95mm,
- > 1 Pc Frosted Glass Stage (For Bottom Illumination)
- > 1 Pc Plastic Black/White Plate (For Top Illumination)
- > Illumination Transmit & Reflected

**Packaging:** Inner Thermocol/Styrofoam  
Packing and Outer Cardboard Box.

