



P-310 Pocket pH Meter

Main Features

- Mini and ultra-thin design with size 180×50×20mm. 10mm diameter electrode can be used to test pH value in tubes.
- Automatic calibration, automatic switch off, low voltage alarm and electrode failure warning.
- IP57 waterproof and easy-to-read LCD screen, replaceable electrode.

Specifications

Model	P-310
Measuring Range	-1.00~15.00 pH
Resolution	0.01pH
Accuracy	±0.2pH
Temperature Compensation	0 ~ 60 ℃
Calibration	1, 2, 3 point
Power	CR2032 button battery
Battery Life	100 hours
IP Rating	IP57
Size & Weight	180×50×20mm/80g(pH Meter) 200×65×50mm/140g(Packing)

Standard Configuration

- | | |
|---|--------|
| 1. P-310 pocket pH meter | 1 unit |
| 2. Backup button battery | 2 pcs |
| 3. pH calibration buffers(pH4.00, pH6.86, pH9.18) | 3 bags |
| 4. Operating manual | 1 copy |

