

- > Ergonomic Design
- > Exceptional Precision
- > Accurate Measurements
- > Highest Industry Standard
- > Outstanding Optical Performance

- > All New Design
- > Binocular Head having Hard coated prisms
- > Quadruple revolving Ball Bearing Nose piece
- > Double layered and Fine Focussing system
- > Special graduated mechanical stage with low positioned coaxial knobs
- > Bright LED Light with Solid state power supply circuitry.
- > Plan Achromatic objectives, Extra Wide Field Eye Piece (Hard Coated) along with high quality condenser lens for better results for all magnifications.

Objectives: 4x, 10x, 40x, 100x (oil)

Eye pieces : EWF 10x (18mm Paired)

Packaging: Inner Thermocol/Styrofoam Packing and Outer Cardboard Box.

Optional Packaging:

- 1) Wooden Box outer cabinet with lock and key (Chargeable).

RESEARCH BINOCULAR MICROSCOPE BM-6bi

ultra



BM-6TR ULTRA

