

GRUN JUNIOR Plus Horizontal Immersed Gel Electrophoresis Unit



Specification:

Make: ORANGE

Cat:GRUN11H

Dimensions: 260x115x75mm

Buffer Volume: 250ml

Number of samples: 11

Small gel tray:60x60mm (2 Nos)

11/11 wells 1mm and 1.5mm thick comb (2 Nos)

8/6 wells 1.5mm thick comb (1 No)

Buffer chamber with safety lid

Positive and negative platinum electrodes

Gel casting stand(77x77x16.5mm): 1 No

Connecting cord: Red and Black

Instruction Manual

Warranty Card