

- > Both Incident & Transmitted lights are provided for ore Microscopy
- > Binocular head replacable with Monocular , and Trinocular head. Bertrand lens, Sliding Analyzer 1 / 2 & 1 / 4 wave compensator.
- > Co-axial Focussing system
- > Rotatable stage divided into 0-360° with vernier
- Adjustable Condensor N.A. 1.25.
- > Polarizer duly graduated with in- out system. Built in light arrangement with 6V-20W halogen lamp.
- > Solid State light control arrangement.
- > Paired eye piece 5x(cross), 10x(ocular), 15x (optional).
- > Objectives P5x, P10x, P40x/45x

Packaging: Inner Thermocol/Styrofoam Packing and Outer Cardboard Box.

ORE MICROSCOPE

OR-01



**OR-02 with
Trinocular Head**