

Power Supply Basic Dual



Specification:

Make : ORANGE

Model : Power Supply Basic Dual

Cat.No : APS2502

Output range : 50,100,150,200,250V/200mA

Type of output : Fixed voltage

Display : LED Indications for Output Fixed Voltages

Output terminals : 2 Pair of recessed banana jacks

Enclosures : Mild Steel with Powder Coated

Operating conditions : 0–40°C, 0–90% humidity in the absence of condensation

Control : Knob Control to set the Voltages

Input protection : 1Amp. Fuse controlled

Input power (actual) : 220 VAC/50hz



APS2502