

SYSTONIC INDIA PRIVATE LIMITED

Fully Automatic Colorimeter S-9121



We are experienced manufacturers and suppliers of Auto Colorimeter for our clients based globally. Manufactured with excellent raw material and technology, they are properly customized and are available in multiple designs. Apart from this, we deliver our products on time.

Technical Specification:

Model	S-9121
Range	400 - 700 nm
Accuracy	± 0.01 O.D.
Stability	± 0.02 O.D. Per Hour
Light Source	Highly Accurate Laser Source
Volume	1 ml. Minimum
Dimension	195X215X95 mm (LXBXH)
Battery Backup	6 V DC Re-chargeable Battery (optional)
Display	LED (7 Segment Display)
Display Model	Absorbance & Wavelength
Power Supply	230V, 10% AC, 50Hz
Filters	Built in 9 Digital Filters 50,000 hours Life
Filter Wavelength	8 Filters- 400, 420, 480, 500, 520, 540, 620, 680 nm
Weight	1.5 Kg. Approx.

Silent Features: -

- Accurate, Direct Readout of Absorbance
- Rugged, Reliable & Easy to Operate
- Long Life Photo Diode Detector
- 1 ml Sample Size

Accessories: -

1. Test Tube (Set of 5)
2. Dust Cover
3. Instruction Manual
4. Adapter