



ultra Series

- > All New Design
- > Monocular Head having hard coated prisms
- > Quadruple revolving Ball Bearing Nose piece
- > Separate Coarse and Fine Focussing system
- > Double layered mechanical stage with low positioned coaxial knobs
- > Bright LED Light with Solid state power supply circuitry.
- > Plan Achromatic objectives, Extra Wide Field Eye Piece (Hard Coated) along with high quality condenser lens for better results for all magnifications.

Objectives : 4x, 10x, 40x, 1 00x (oil)

Eye pieces: EWF 10x (18mm Paired)

Packaging: Inner Thermocol/Styrofoam Packing and Outer Cardboard Box.

MONOCULAR MICROSCOPE BM-6mo

ultra

