

UV TRANSILLUMINATOR



Make : ORANGE

Model : UV Transilluminator BLB Research

Cat.No : UVG3652

UV source :: 365nm wave length

Viewing area : 25x25cm

UV Tube : 10x8 watts UV Tube

Housing : Anti Corrosion treated with Powder coated Mild Steel

UV Safety Shield : Clear Transparent lid for excellent clarity

Input : 220V AC, 50-60HZ

Fuse : 2.5Amp

Anti Slip rubber foot provided for base

Instruction Manual & Warranty card

RESEARCH

UVG3652

