



Body

Monocular fixed straight tube (Length - 160mm)

Illumination

LED Lamp

Stage

110 x 110mm square stage is provided with stage clips to hold the slide (optional Mechanical Stage).

Focussing

Focussing control knob on both sides of the instrument is conveniently located for easy operation of the rack and pinion with 15mm range.

Nose Piece

Triple revolving Nosepiece with positive click stops

Screen

150mm diameter glass projection screen is provided which can be rotated over 360 degrees on Microscope tube.

Objectives

Achromatic 4x, 10x & 40x

Packaging: Inner Thermocol/Styrofoam Packing and Outer Cardboard Box.

STUDENT
PROJECTION
MICROSCOPE
(Doom Type)

KW-700C

